











Directions for Using the Marvel Violet Ray

1. Insert glass electrode into the instrument, with straight direct push.
2. Attach the connection plug to the lamp socket.
3. Turn the current on.
4. Adjustment of the Violet Ray discharge is obtained by the adjusting knob. To the right the discharge increases, to the left it decreases.
5. Start any treatment with a mild discharge and increase it gradually to the required strength.
6. Standard winding for the "Marvel" is 110 Volts in the U. S. A. (Twenty watts). For lower and higher voltage special windings are necessary. Information will gladly be furnished upon request.
7. Take good care of the connecting cord. Keep it dry and free from knots and twists. Have it replaced in time if signs of wear appear.

Don'ts

1. Do not attempt to repair the apparatus. If the necessity arises, write to the "Eastern Laboratories, Inc., N. Y. C." advising of your difficulties and you will be informed what to do.
2. Do not use the "Marvel" in connection with any tonic or medicament containing a large percentage of alcohol on the hair or on any part of the body.
3. When not in use, machines should always be disconnected from lamp socket.



Note: Damages arising from carelessness are not covered by our guarantee. If instructions are properly followed the "Marvel" will last indefinitely.

[illegible]

No. 8-1—Standard surface electrode for general use \$1.00

No. 8-4—Metal electrode or Batometer, used in giving indirect treatments \$1.25

No. 8-5—Distal throat electrode, also excellent for use on neck, arm, or any curved surface of the body \$1.00

No. 8-3—Large comb electrode for scalp \$1.50

No. 8-2—Large condenser electrode, equipped with lead-in device which carries current right to surface of electrode. Indicated for use in cases of Rheumatism, Lumbago, etc. \$1.00

No. 8-6—Sphulac electrode used for stimulating the entire body by applying to the spine \$1.00

No. 8-20—Pulsarization electrode for removal of warts, moles, discolored growths, etc. \$2.00

No. 8-10—Insulated Ear Electrode \$2.00

No. 8-18—Single eye electrode \$1.00

No. 8-19—Artzels electrode. Develops maximum X-rays to cause chemical pyro effect. Used for treatment of Acan, Vulgaris, Ringworm, Lice, etc. \$1.00

No. 8-16—Large Kathodeion (cathode) electrode. Medicated cotton can be inserted into hollow surface and the high frequency current passes through medicated cotton. Small particles of the medication are thus driven into the tissue by molecular bombardment. The effect is powerful. \$1.00

No. 8-17—Double Eye electrode \$2.00

No. 8-14—Plate rectal electrode, small dia. (1 1/2") \$1.00

No. 8-15—Plate rectal electrode, medium diam. (1 3/4") \$1.00

No. 8-13—Plate rectal electrode, large dia. (1 1/2") \$1.00

No. 8-20—Insulated vaginal electrode \$2.00

No. 8-12—Insulated rectal electrode, small (1 1/2") \$2.00

No. 8-14—Insulated rectal electrode, medium (1 3/4") \$2.00

No. 8-13—Insulated rectal electrode, large (1 1/2") \$2.00

No. 8-21—Insulated vaginal electrode for the generation of ozone within the vagina \$1.00

No. 8-21—Plain vacuum electrode, good for many uses as for instance for treatment of the vocal cords, overworn, or hands for indirect applications. \$1.00

No. 8-1—Insulated Internal throat electrode \$2.00

No. 8-10—Plain Internal throat electrode \$1.00

No. 8-50—Insulated urethral electrode \$1.00

No. 8-51—Plain urethral electrode \$1.00

No. 8-9—Insulated prostatic electrode \$2.00

No. 8-7—Insulated spatulum or tongue electrode \$2.75

No. 8-8—Plain nasal electrode \$1.00

No. 8-11—Insulated nasal electrode \$2.00

No. 8-18—Hollow surface shaped electrode, gentler action within cavity and for work over certain parts, like cheekbones \$1.25

No. 8-19—Core electrode \$1.00

No. 8-12—Kryof electrode \$1.25

No. 8-42—Hollow kataphoretic (facial) electrode, to be used the same as No. 8-14 \$2.50

No. 8-43—Dental electrode for gum treatment \$2.00

No. 8-44—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-45—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-46—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-47—Dental electrode, straight form \$2.00

No. 8-48—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-49—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-50—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-51—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-52—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-53—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-54—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-55—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-56—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-57—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-58—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-59—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-60—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-61—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-62—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-63—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-64—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-65—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-66—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-67—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-68—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-69—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-70—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-71—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-72—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-73—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-74—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-75—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-76—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-77—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-78—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-79—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-80—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-81—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-82—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-83—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-84—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-85—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-86—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-87—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-88—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-89—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-90—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-91—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-92—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-93—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-94—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-95—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-96—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-97—Dental electrode for gum treatment, right side of jaw \$2.00

No. 8-98—Dental electrode for gum treatment, left side of jaw \$2.00

No. 8-99—Dental electrode for gum treatment, inside of jaw \$2.00

No. 8-100—Dental electrode for gum treatment, right side of jaw \$2.00







